



Wheelchair Power Attachment



2025
CATALOG



Email: info@elesmart.com

Tel: +86-13685249056

Web: www.elesmart.com

Add: West Industrial Park, Hutang, Wujin District, Changzhou City, Jiangsu Province, China

ELES MART CO.,LTD

ELESMART-A4 **03**

01 ELESMART-A3

05 ELESMART-A5

ELESMART-A6 **07**

09 ELESMART-A8

ELESMART CO.,LTD



HOT SELLING COUNTRIES



CE CERTIFICATION












GOOD QUALITY

About us

ELESMART specializes in the development and production of wheelchair power attachment ,electric handbikes. Founded in 2011, the company is supported by a team of six technical experts, each with over 25 years of experience in the medical equipment industry. Our mission is to continuously enhance mobility for individuals with disabilities, enabling them to integrate into society and achieve a better quality of life.

At ELESMART, improving the daily mobility and quality of life for wheelchair users through our products is our ultimate motivation and driving force.



 Wheel Diamter	12"	 Torque	40Nm	 Display	COLOR TFT+USB
 Motor Power	36V 350W	 Lithium Battery	36V 7.8/12.8Ah	 Frame	Standard
 Speed*	25KM/H	 Range (KM)*	40/60	 Connection Method	Side


- 1.Gear and Speed Parameters: 3 forward gears, maximum speed of 25KM/H, 1 reverse gear, right-hand control for forward motion.
- 2.Cruise Control Function: Enabled by default, can be turned off.
- 3.Maximum Speed: Affected by actual user conditions, load, road conditions, etc., with a variation of ± 5 KM.
- 4.Range: Measured with a weight of 50kg work on 1 level on a flat road. Range depends on user weight, road conditions, wheelchair tire condition, battery condition, speed, and driving habits.
- 5. Charger: Complies with EU regulations, EMC tested, 36V dimensions is 104.3*53*26mm, 300g, Easy to take !

Colours



Wheel Diamter	10"	Torque	40Nm	Display	LCD4
Motor Power	36V 350W	Lithium Battery	36V 7.8Ah	Frame	Foldable
Speed*	25KM/H	Range (KM)*	40	Connection Method	Side




- 1. Gear and Speed Parameters: 3 forward gears, maximum speed of 25KM/H, 1 reverse gear, right-hand control for forward motion.
- 2. Cruise Control Function: Enabled by default, can be turned off.
- 3. Maximum Speed: Affected by actual user conditions, load, road conditions, etc., with a variation of ± 5 KM.
- 4. Range: Measured with a weight of 50kg work on 1 level on a flat road. Range depends on user weight, road conditions, wheelchair tire condition, battery condition, speed, and driving habits.
- 5. Charger: Complies with EU regulations, EMC tested, 36V dimensions is 104.3*53*26mm, 300g, Easy to take !



Colours



	Wheel Diamter	12"
	Motor Power	36V/48V 500W
	Speed*	25KM/H

	Torque	45Nm
	Lithium Battery	36V 10.4/12.8/16Ah
	Range (KM)*	50/60/70

	Display	COLOR TFT+USB
	Frame	Suspension
	Connection Method	Side

- 1. Gear and Speed Parameters: 3 forward gears, maximum speed of 25KM/H, 1 reverse gear, right-hand control for forward motion.
- 2. Cruise Control Function: Enabled by default, can be turned off.
- 3. Maximum Speed: Affected by actual user conditions, load, road conditions, etc., with a variation of ± 5 KM.
- 4. Range: Measured with a weight of 50kg work on 1 level on a flat road. Range depends on user weight, road conditions, wheelchair tire condition, battery condition, speed, and driving habits.
- 5. Charger: Complies with EU regulations, EMC tested, 36V dimensions is 104.3*53*26mm, 300g, Easy to take !



Colours



	Wheel Diameter	12"
	Motor Power	36V/48V 500W
	Speed*	25KM/H

	Torque	45Nm
	Lithium Battery	36V 10.4/12.8/16Ah
	Range (KM)*	50/60/70

	Display	COLOR TFT+USB
	Frame	Integrated
	Connection Method	Center

- 1. Gear and Speed Parameters: 3 forward gears, maximum speed of 25KM/H, 1 reverse gear, right-hand control for forward motion.
- 2. Cruise Control Function: Enabled by default, can be turned off.
- 3. Maximum Speed: Affected by actual user conditions, load, road conditions, etc., with a variation of ± 5 KM.
- 4. Range: Measured with a weight of 50kg work on 1 level on a flat road. Range depends on user weight, road conditions, wheelchair tire condition, battery condition, speed, and driving habits.
- 5. Charger: Complies with EU regulations, EMC tested, 36V dimensions is 104.3*53*26mm, 300g, Easy to take !



Colours



	Wheel Diameter	16"
	Motor Power	48V 1000W
	Speed*	40KM/H

	Torque	85Nm
	Lithium Battery	48V 15/20Ah
	Range (KM)*	80/100

	Display	COLOR TFT+USB
	Frame	Integrated
	Connection Method	Center

- 1. Gear and Speed Parameters: 3 forward gears, maximum speed of 40KM/H, 1 reverse gear, right-hand control for forward motion.
- 2. Cruise Control Function: Enabled by default, can be turned off.
- 3. Maximum Speed: Affected by actual user conditions, load, road conditions, etc., with a variation of ± 5 KM.
- 4. Range: Measured with a weight of 50kg work on 1 level on a flat road. Range depends on user weight, road conditions, wheelchair tire condition, battery condition, speed, and driving habits.
- 5. Charger: Complies with EU regulations, EMC tested, 36V dimensions is 104.3*53*26mm, 300g, Easy to take !